

# STERIWAY

## Amalgam Separator




Equipped and Committed to assist the dental industry achieve environmentally sound and regulatory compliant practices.

provides the complete solution:

- ISO 11143 certified at greater than 99%
- Professionally engineered, robust, simple design.
- Easy to operate and maintenance free.
- Flexible design for hard-to-fit applications.
- Functional for wet or dry vacuum systems (install before pump on wet vacuum system, install before tank on dry vacuum system).

Simple container change process includes:

- NO TOOLS, NO MESS
- PACKAGING FOR DELIVERY TO CERTIFIED RECYCLER
- ONLINE CERTIFICATE ISSUED FOR PRACTICE



dta-med.com

دانش طب امروز  
۸۸۹۳۰۹۲۷



[www.dta-med.com](http://www.dta-med.com)

[info@dta-med.com](mailto:info@dta-med.com)

